		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report	
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989
Originating Department Analytics	Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:12:59 PM UTC • Page 1 of 23

Validation of Surveillance Assay for SARS-CoV-2 Testing of Pooled Patient Swab Samples

(The Infinity BiologiX TaqMan SARS-CoV-2, Flu A, Flu B RT-PCR Surveillance Assay will be performed in the Infinity BiologiX LLC Clinical Genomics Laboratory, a Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, certified high-complexity laboratory).

Test Summary:

The Infinity BiologiX TaqMan SARS-CoV-2, Flu A, Flu B RT-PCR Assay is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test intended for the qualitative detection of nucleic acid from SARS-CoV-2, Influenza A, and Influenza B in saliva, oropharyngeal (throat) swab, nasopharyngeal swab, anterior nasal swab, mid-turbinate nasal swab, and bronchoalveolar lavage (BAL) fluid from individuals suspected of COVID-19 and/or Influenza A/B

Testing is limited to Infinity BiologiX LLC locations that are certified under Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, and meets requirements to perform high-complexity tests.

Results are for the detection and identification of SARS-CoV-2, Influenza A, and Influenza B RNA for surveillance purposes. The SARS-CoV-2, Influenza A, or Influenza B RNA are generally detectable in respiratory specimens collected using nasopharyngeal swabs during the acute phase of infection.

Specimen pooling for surveillance purposes can increase throughput and conserve testing resources, but validation is needed to ensure that the reduced sensitivity due to pooling does not increase the false-negative rate. We use of ThermoFisher's swab samples in a pooled fashion. Pooled specimen testing can help increase testing capacity for SARS-CoV-2 surveillance with a low risk of false-negative results. In a simple pooled testing scheme, a number of individually collected specimens are combined in a single well or tube and tested together. If the pooled test result is negative, results for all individual specimens may be immediately reported as negative. If a pool is positive, then each individual in the pool will be asked to quarantine until the individual testing results are available. The optimal number of specimens that can be included in a pool to maximize efficiency is determined by the prevalence of positive specimens in the population being tested and is further constrained by the sensitivity of the test to reliably detect a positive signal in a diluted negative specimen pool. It is critical, therefore, to validate pooling strategies for diagnostics tests to ensure that the false-negative rate remains below an acceptable threshold

Testing with the TaqMan SARS-CoV-2/Flu A, Flu B RT-PCR test is intended for use by qualified laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures. The assay is intended for testing pooled patient samples for surveillance purposes.

Assay Protocol Overview:


The Infinity BiologiX TaqMan SARS-CoV-2, Flu A, Flu B, RT-PCR Assay is a real-time reverse transcription assay designed to detect RNA from SARS-CoV-2, Flu A, Flu B in nasopharyngeal swabs collected from participants. The purpose of this validation is to enable testing of pooled COVID-19/Flu A/Flu B nasopharyngeal swab samples for surveillance purposes

Swab specimens must be collected, transported and stored using the Miraclean foam tip clean room swabs (ThomasScientific, Catalog # **20A00M411 / MFS-741**). Swab specimens must be transported dry and stored at ambient temperature and tested within 72 hours of collection when stored at ambient temperature.

Swab Pooling, Virus spike-in, and Reconstitution strategy

Swab Pooling and virus spike-in

Twenty-four (24) swab samples were collected using manufacturer's instructions from twenty-four (24) patients tested negative for SARS-CoV-2, which was confirmed via IBX SARS COV-2 RT-PCR assay. All 24 swabs collected from patients

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:12:59 PM UTC • Page 2 of 23

are placed in one 50 ml Falcon tube. A 25th swab is then spiked with COVID19 virus (*BEI Resources* NR-52287, Lot # 70035888) at the concentrations representing pools with low viral concentrations (1,500 copies/ml – 3x Limit of Detection (LoD)), or weak positive sample in a pool (1,500 copies/mL, or strong positive sample in a pool (1,312,500 copies/mL, diluted in 5–10 µl phosphate buffered saline (PBS) and then the spiked swab is added to the same tube containing 24 negative swabs. Different swab pools were created to determine the performance across a wide variety of pooled specimens as laid out in **Table 1** below.

Table 1: Strategy for creating pooled swab samples

Performance Criteria	Transcript copies/mL	Number of negative swabs per pool	Number of positive swabs per pool	Total number of pools	Total number of swabs
Confirmation of performance in a negative pool	0	25	0	9	225
Confirmation of performance at low viral concentrations	1,500 (without clinical matrix)	24	1	9	225
Detection of weak positive sample in a pool	1,500 (with clinical matrix)	24	1	12	300
Detection of strong positive sample in a pool	1,312,500	24	1	6	150
	Total			36	900

Reconstitution of dry swab samples:

Prior to testing the pooled specimens using TaqMan SARS-CoV-2, Flu A, Flu B RT-PCR Assay, the dry pooled swabs are reconstituted in 6 ml of PBS followed by vortexing for 30 s with intermittent pulsing. The reconstituted samples are then incubated for 20 min at room temperature. After incubation, the samples are vortexed again for 30 s with intermittent pulsing. 300 µl of the reconstituted sample is used in subsequent steps.


RNA extraction for all swab samples is performed using the PerkinElmer Chemagic 360 automated specimen processing system with the Chemagic Viral DNA/RNA 300 Kit H96. The input sample volume is 300µL, the elution volume is 60µL. Reverse transcriptase-PCR (RT-PCR) is performed using the Applied Biosystems TaqMan SAS-CoV-2/Flu A/Flu B assay with 12.5µL of the extracted sample.

Instruments used for the test:

The Infinity BiologiX TaqMan SARS-CoV-2, FluA, FluB Assay is for use with the ThermoFisher Applied Biosystems QuantStudio 5 Real-Time PCR System equipped with software v1.3, or the Applied Biosystems ViiA7 Real-Time PCR System with the Applied Biosystems QuantStudio 5 software v1.3 for data analysis, and Perkin Elmer Chemagic 360 extraction instrument (software v6.3.0.3).

Table 2: Reagents and materials required for use of the Infinity BiologiX TaqMan SARS-CoV-2, FluA, FluB Assay

Reagents	Manufacturer	Catalog #
Chemagic Viral DNA/RNA 300 Kit H96	PerkinElmer	CMG-1033-S
96 well Deep Well Plates	PerkinElmer	43001-0120
TaqMan SARS-CoV-2/Flu A/Flu B Assay	ThermoFisher Scientific	A47701
384 well PCR plate	ThermoFisher Scientific	4483273
Optical adhesive PCR plate cover	ThermoFisher Scientific	4311971
Nuclease-free water	--	--
Ethanol (96-100%)	--	--

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report	
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989
Originating Department Analytics	Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:12:59 PM UTC • Page 3 of 23

Controls:

The controls supplied with the ThermoFisher TaqMan COVID-19 Combo Kit are described in **Table 3**.

Table 3. Controls supplied with the Applied Biosystems TaqMan COVID-19 Combo Kit

Control Type	Purpose	Frequency of Testing
Negative	To monitor for cross-contamination during RNA extraction and RT-PCR	Once per batch of specimens
Positive	To monitor the integrity of the RT-PCR reagents and process	Once per run of RT-PCR
Internal (MS2 Phage)	To monitor the integrity of nucleic acid extraction and RT-PCR for each specimen	Added to each specimen and the Negative Control prior to extraction

In addition to these controls, a No Template Control containing none of the SARS-CoV-2 targets or the Internal Control is included in every PCR run. The results from the controls are interpreted according to the criteria shown in **Table 4**. If the results obtained with the Positive, Negative and No Template Controls do not meet the criteria shown, the results from the entire batch of samples are considered invalid and repeat testing must be performed.

Table 4: Ct values for controls that must be observed to obtain valid results


Control	Ct value (Optical Channel)	
	N/S Gene	MS2 Phage
	(VIC)	(JUN)
Negative	>40	<37
Positive	<37	Undetermined
No Template	Undetermined	Undetermined
Internal	Any	Any

Interpretation of Results:

The results from testing of pooled swab samples are interpreted according to the criteria described in **Table 5**.

Table 5: Result interpretation of pooled patient samples

Ct Value (Optical Channel)		Result Interpretation
N/S Gene	MS2 Phage	
(VIC)	(JUN)	
Undetermined	<37	Negative
<37	<37	Positive
<37	Undetermined	Re-test
Undetermined	Undetermined	Re-test

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report		
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989	
Originating Department Analytics		Affected Department(s) Analytics	Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:12:59 PM UTC • Page 4 of 23

Reportable range

The TaqMan SARS-CoV-2/Flu A/Flu B Assay for detection of SARS-CoV-2/FluA/FluB in pooled clinical swab samples is a qualitative test and therefore a reportable range does not apply.

Reference range

The reference range will be "SARS-CoV-2/FluA/FluB not detected" and will not be validated by a clinical study.

Performance Evaluation:

Accuracy

To assess the test accuracy of detecting negative and positive contrived pooled swab specimens, we ran the TaqMan SARS-CoV-2, Flu A, Flu B RT-PCR Assay on four different swab pools as laid out in the **Table 6** below. After reconstituting the 25 dry swab pools in 6 mL of PBS, we obtained 10 extraction replicates of 300 µL each from the pool and ran RT-PCR on each replicate. 90/90 (100%, 95% CI [95.98% - 100%]) negative pooled swab specimens tested negative. 90/90 (100%, 95% CI [95.98% - 100%]) pooled swab specimens 97% - 100%) pooled swab specimens with clinical matrix spiked with 1,500 copies/mL of the virus tested positive and 59/60 (95%, 95% CI [91% - 99.96%]) pooled swab specimens with clinical matrix spiked with 1,312,500 copies/mL of the virus tested positive. The results of the assay are summarized in **Table 6** below. The positive (PPA) and negative percent agreement (NPA) are summarized in **Table 7** below.


Table 6: Summary results of pooled swabs. *One pooled sample from the strong positive pooled swabs group tested inconclusive.

Performance Metric	Transcript Copies/mL	Pool type	Swab Pools	Analysis	Positive	Negative
Negative pool	0	Negative	9	Negative (%)	0 (0)	90 (100)
Confirmation of performance at low viral concentration	1,500	Positive	9	Positive (%)	90 (100)	0 (0)
Detection of weak positive swabs in a pool	1,500	Positive	12	Positive (%)	120 (100)	0 (0)
Detection of strong positive swabs in a pool	1,312,500	Positive	6	Positive (%)	59 (98.33)	1* (1.67)

Table 7: Positive percent (PPA) and negative percent agreement (NPA) for Pooled swab samples vs expected results tested across three different days with three sets of independent samples

Samples Tested in 25-swab pool	Expected Results	
	Positive	Negative
Negative	0 (0%)	90 (100%)
1,500 (without clinical matrix)	90 (100%)	0 (0)
1,500 (with clinical matrix)	120 (100%)	0 (0)
1,312,500 (with clinical matrix)	59 (98.33%)	1 (1.67%)
All Positive	269 (99.62%)	1 (0.37%)

Acceptance criteria: > 95% percent agreement between replicate samples. **PASS**

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:12:59 PM UTC • Page 5 of 23

Precision

Repeatability

To assess the repeatability (within-run precision) of pooled swab specimen testing, we performed ten replicate measurements of the contrived validation samples described in the section on Accuracy within a single batch. We then assessed the concordance between the ten replicate measurements and determined the percent of concordant results.

In total, 359/360 (98.33%, 95% CI [98.46% - 99.99%]) replicates were concordant between ten repeat measurements run on the same batch for 36 swab pools. See Table 8, 9, and 10 for between run variations across 3 different days. See **Appendix** for detailed run data.


Acceptance criteria: > 95% percent agreement between duplicates. **PASS**

Table 8: Repeatability of pooled swab specimen testing (within in run precision)

	Performance Criteria	Copies/mL	n	Analysis	Day 1		Day 2		Day 3	
					Target (Optical Channel)		Target (Optical Channel)		Target (Optical Channel)	
					N/S Gene	MS2	N/S Gene	MS2	N/S Gene	MS2
					(VIC)	(JUN)	(VIC)	(JUN)	(VIC)	(JUN)
25-Swab Pools	Performance in a negative pool	0	30	Mean	NA	19.01	NA	21.97	NA	22.84
				SD	NA	1.15	NA	0.31	NA	0.67
				CV	NA	6%	NA	1%	NA	3%
	Performance at low viral concentration	1,500	30	Mean	31.66	18.94	33.56	22.16	32.62	22.58
				SD	2.94	1.19	0.61	0.22	2.40	0.42
				CV	9%	6%	2%	1%	7%	2%
	Detection of weak positive swabs in a pool	1,500	40	Mean	28.83	19.57	30.93	22.19	31.96	22.70
				SD	2.88	0.81	0.93	0.33	2.53	0.52
				CV	10%	4%	3%	1%	8%	2%
Detection of strong positive swabs in a pool	1,312,500	20	Mean	20.69	20.37	20.90	22.22	22.29	22.89	
			SD	1.54	0.95	0.27	0.36	1.25	0.68	
			CV	0.07	0.05	0.01	2%	6%	3%	

Table 9: Day 1 Positive percent (PPA) and negative percent agreement (NPA) for Pooled swab samples vs expected results

Samples Tested in 25-swab pool	Expected Results	
	Positive	Negative
Negative	0 (0%)	30 (100%)
1,500 (without clinical matrix)	30 (100%)	0 (0)
1,500 (with clinical matrix)	40 (100%)	0 (0)
1,312,500 (with clinical matrix)	19 (95.00%)	1 (5%)
All Positive	89 (98.88%)	1 (1.12%)

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report		
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989	
Originating Department Analytics		Affected Department(s) Analytics	Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 6 of 23

Table 10: Day 1 Performance metrics

Performance metric	Description	Performance
Confirmation of performance at low viral concentration	Percent of swab pools positive at 1500 copies/mL (without clinical matrix)*	30 / 30 = 100%; 95% CI [88.4%-100%]
Detection of weak positive swabs in a pool	Percent of swab pools positive at 1500 copies/mL (with clinical matrix)*	40 / 40 = 100%; 95% CI [91.2%-100%]
Detection of strong positive swabs in a pool	Percent of swab pools positive at 1,312,500 copies/mL (with clinical matrix)*	19 / 20 = 95%; 95% CI [75.1%-100%]

Table 11: Day 2 Positive percent (PPA) and negative percent agreement (NPA) for Pooled swab samples vs expected results


Samples Tested in 25-swab pool Transcript Copies/mL	Expected Results	
	Positive	Negative
Negative	0 (0%)	90 (100%)
1,500 (without clinical matrix)	30 (100%)	0 (0)
1,500 (with clinical matrix)	40 (100%)	0 (0)
1,312,500 (with clinical matrix)	20 (100%)	0 (0)
All Positive	90 (100%)	0 (0)

Table 12: Day 2 Performance metrics

Performance metric	Description	Performance
Confirmation of performance at low viral concentration	Percent of swab pools positive at 1500 copies/mL (without clinical matrix)*	30 / 30 = 100%; 95% CI [88.4%-100%]
Detection of weak positive swabs in a pool	Percent of swab pools positive at 1500 copies/mL (with clinical matrix)*	40 / 40 = 100%; 95% CI [91.2%-100%]
Detection of strong positive swabs in a pool	Percent of swab pools positive at 1,312,500 copies/mL (with clinical matrix)*	20 / 20 = 100%; 95% CI [83.2%-100%]

Table 13: Day 3 Positive percent (PPA) and negative percent agreement (NPA) for Pooled swab samples vs expected results

Samples Tested in 25-swab pool Transcript Copies/mL	Expected Results	
	Positive	Negative
Negative	0 (0%)	90 (100%)
1,500 (without clinical matrix)	30 (100%)	0 (0)
1,500 (with clinical matrix)	40 (100%)	0 (0)
1,312,500 (with clinical matrix)	20 (100%)	0 (0)
All Positive	90 (100%)	0 (0)

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report		
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989	
Originating Department Analytics		Affected Department(s) Analytics	Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 7 of 23

Table 14: Day 3 Performance metrics

Performance metric	Description	Performance
Confirmation of performance at low viral concentration	Percent of swab pools positive at 1500 copies/mL (without clinical matrix)*	30 / 30 = 100%; 95% CI [88.4%-100%]
Detection of weak positive swabs in a pool	Percent of swab pools positive at 1500 copies/mL (with clinical matrix)*	40 / 40 = 100%; 95% CI [91.2%-100%]
Detection of strong positive swabs in a pool	Percent of swab pools positive at 923,000 copies/mL (with clinical matrix)*	20 / 20 = 100%; 95% CI [83.2%-100%]

Reproducibility

To assess the reproducibility (between-run precision) of pooled saliva specimens, we performed three measurements of the clinical validation samples described in Accuracy on three separate days performed by three different operators on three different instruments. We then assessed the concordance between the three runs and determined the percent of concordant results to determine between run variability


In total, 119/120 (99.17%; 95% CI [95.44% - 99.98%]) replicates run each day were concordant between ten measurements run within the same batch across three different days. See Table 8, 9, and 10 for between run variations across 3 different days. See **Table 15** for summary of the results for assay reproducibility and **Appendix** for detailed run data.

Acceptance criteria: > 95% concordance of samples between runs. **PASS**

Table 15: Reproducibility of pooled swab specimen testing (between run precision)

	Performance Criteria	Description	Copies/mL	n	Analysis	Target (Optical Channel)	
						N/S Gene	MS2
						(VIC)	(JUN)
25-Swab Pools	Confirmation of performance in a negative pool	Pooled sample without virus spike-in representing a negative pool	0	90	Mean	NA	21.27
					SD	NA	1.82
					CV	NA	9%
	Confirmation of performance at low viral concentration	Percent of swab pools positive at 1,500 copies/mL (without clinical matrix)*	1,500	90	Mean	32.61	21.23
					SD	2.33	1.79
					CV	7%	8%
	Detection of weak positive swabs in a pool	Percent of swab pools positive at 1,500 copies/mL (with clinical matrix)*	1,500	120	Mean	30.57	21.49
					SD	2.61	1.50
					CV	9%	7%
	Detection of strong positive swabs in a pool	Percent of swab pools positive at 1,312,500 copies/mL (with clinical matrix)*	1,312,500	60	Mean	21.31	21.85
					SD	1.34	1.27
					CV	6%	6%

Analytical Sensitivity

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report		
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989	
Originating Department Analytics		Affected Department(s) Analytics	Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 8 of 23


The Limit of Detection (LoD) is defined as the lowest SARS-CoV-2, influenza A, influenza B RNA viral concentration that is successfully detected with a probability of 95% or greater. The LoD was determined using *in vitro* heat-inactivated SARS-CoV-2 virus (BEI NR-52287), Flu A (Zeptomatrix, A/Perth/16/2009 and A/Brisbane/59/2007), Flu B (Zeptomatrix, B/Florida/04/2006, and B/Wisconsin/01/2010) that was diluted in SARS-CoV-2, FluA, FluB negative single nasopharyngeal swab matrix. An initial estimate of the LoD with the Applied Biosystems QuantStudio 5 Real-Time PCR System was obtained by testing three replicates at each of four different target levels: 1000, 500, 200 and 100 copies/mL. The lowest level at which all three replicates were positive for SARS-CoV-2, FluA, Flu B targets was determined. The estimated LoD was confirmed by testing an additional 20 replicates at the same target level. All 20 NP swab replicates produced the expected results. The LoD was therefore confirmed to be as that shown in **Table 16** below.

Table 16: Limit of Detection of Infinity BiologiX TaqMan SARS Cov-2, FluA, FluB, combo assay

Target	Strain/details	Vendor	Catalog #	Lot Number	Stock Viral GCE/mL	LoD in GCE/mL
SARS-CoV-2	Isolate USA-WA1/2020	BEI	NR-52287	70035888	1.75E+09	100
Flu A	A/Perth/16/2009	ZeptoMetrix	0810251CF	324079	3.11E+09	200
Flu A	A/Brisbane/59/2007	ZeptoMetrix	0810244CF	323919	1.58E+11	500
Flu B	B/Florida/04/2006	ZeptoMetrix	0810255CF	312479	1.07E+10	500
Flu B	B/Wisconsin/01/2010	ZeptoMetrix	0810241CF	323015	1.37E+10	1000

Furthermore, an initial estimate to determine the limit of detection in a pooled swab sample with the Applied Biosystems QuantStudio 5 Real-Time PCR System was obtained by testing three replicates of 25 dry pooled swab samples (24 negative samples, and one sample spiked at 1,200 or 1,500 copies/mL). The lowest level at which all three pooled replicates were positive for SARS-CoV-2 N/S Gene was 1,500 copies/mL. All 20 replicates produced the expected results for each SARS-CoV-2 target, and the LoD was therefore confirmed to be 1,500 copies/mL for 25 dry swabs per pool specimen.

Validated results of the 25 dry swabs per pool across three days are summarized in **Table 17** below.

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report		
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989	
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 9 of 23

Table 17. Validation of pooling test performance with 25 dry swabs per pool across three days

Performance Criteria	Description	Copies /mL	Swab pools	Analysis	Target (Optical Channel)	
					N/S Gene	MS2
					(VIC)	(JUN)
25 - Swab Pool	Confirmation of performance in a negative pool	0	9	Negative (%)	0 (0)	90 (100)
				Mean (SD)	NA (NA)	21.27 (1.82)
				CI (95%)	96% - 100%	
	Confirmation of performance at low viral concentration	1,500	9	Positive (%)	90 (100)	90 (100)
				Mean (SD)	32.61 (2.33)	21.23 (1.79)
				CI (95%)	96% - 100%	
	Detection of weak positive swabs in a pool	1,500	12	Positive (%)	120 (100)	120 (100)
				Mean (SD)	30.57 (2.60)	21.48 (1.49)
				CI (95%)	97% - 100%	
	Detection of strong positive swabs in a pool	1,312,500	6	Positive (%)	59 (98.33)	1 (1.67)
				Mean (SD)	21.30 (1.33)	21.84 (1.27)
				CI (95%)	91.06% - 99.96%	

*Virus used to prepare swab samples: BEI Resources NR-52287, Lot # 70035888

Analytical Specificity


Inclusivity

The Infinity BiologiX TaqPath SARS-CoV-2 Assay is a modification of the previously authorized ThermoFisher Applied Biosystems TaqPath COVID-19 Combo Kit. The assay targets specific genomic regions of the SARS-CoV-2 nucleocapsid (N) gene, spike (S) gene, and ORF1ab region, influenza A and Influenza B. Inclusivity was demonstrated by mapping the primers and probes to complete SARS-CoV-2, Influenza A, Influenza B genomes that were available in the GenBank and GISAID (Global Initiative on Sharing All Influenza Data) databases as of March 5, 2020. For all primers and probes, there was 100% homology to each of the SARS-CoV-2, Flu A, and Flu B sequences analyzed

Acceptance criteria: Greater than 80% homology with each sequence analyzed. **PASS**

Cross reactivity

The Infinity BiologiX TaqMan SARS-CoV-2/Flu A/Flu B Assay is a modification of the previously authorized ThermoFisher Applied Biosystems TaqMan COVID-19 Combo Kit. In silico analysis was performed to predict cross-reactivity between the TaqPath™ RT-PCR COVID-19, Flu A, Flu B Combo kit (COVID/Flu test) primer/probe sequences and a series of microbes commonly found in respiratory specimens per [the] TaqPath™ RT-PCR COVID, Flu A, Flu B Kit Cross-Reactivity Study Protocol. Based upon BLAST analysis, 33 isolates representing 21 total bacterial and fungal organisms as shown in **Table 18** had >=80% homology to either the reverse primer, forward primer or probe sequence of the COVID/Flu test SARS-


		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report	
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989
Originating Department Analytics	Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 10 of 23

CoV-2 N gene, Influenza A, or MS2 internal positive control assays; no primer or probe homology of $\geq 80\%$ was identified with the Influenza B assay. For the select isolates that shared $\geq 80\%$ identity with a given assay primer or probe sequence, sequence identity with other assay components was not detected. Given that both assay primers are required to produce an amplicon and all three assay components are required to generate signal, it was determined that there is little to no chance of cross-reactivity between assays and other respiratory microbes, or interference with detection of SARS-CoV-2, influenza A and influenza B viruses, or the MS2 internal control. Based on this analysis, significant amplification of non-target sequences that could result in cross-reaction (false-positive results) or interference (false-negative results) was considered unlikely to occur.

Because cross reactivity was assessed for the SARS-CoV-2 S gene assay used in this test under [under previous EUA for TaqPath COVID-19 assay], it was not required to evaluate the SARS-CoV-2 S gene assay in this study.

Acceptance criteria: 0% Cross reactivity or interference with detection. **PASS**

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report		
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989	
Originating Department Analytics		Affected Department(s) Analytics	Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 11 of 23


Table 18: Organisms and viruses evaluated for potential cross-reaction and/or interference with the TaqMan-SARS-Cov-2/Flu A/Flu B RT-PCR assay

<i>Bacillus anthracis</i>	<i>Adenovirus</i>
<i>Bordetella pertussis</i>	<i>Enterovirus</i>
<i>Chlamydia psittaci</i> ^[1]	<i>Human Coronavirus 229E</i>
<i>Chlamydomphila pneumoniae</i>	<i>Human Coronavirus HKU1</i>
<i>Corynebacterium diphtheriae</i>	<i>Human Coronavirus NL63</i>
<i>Coxiella burnetii</i> ^[1]	<i>Human Coronavirus OC43</i>
<i>Haemophilus influenzae</i>	<i>Human Metapneumovirus</i>
<i>Legionella non-pneumophila (Legionella longbeachae)</i>	<i>Influenza C virus</i>
<i>Legionella pneumophila</i>	<i>MERS Coronavirus</i>
<i>Leptospira spp.</i>	<i>Parainfluenza 1 virus</i>
<i>Moraxella catarrhalis</i>	<i>Parainfluenza 2 virus</i>
<i>Mycobacterium tuberculosis</i>	<i>Parainfluenza 3 virus</i>
<i>Mycoplasma pneumoniae</i>	<i>Parainfluenza 4 virus</i>
<i>Neisseria elongata</i>	<i>Parechovirus</i>
<i>Neisseria meningitidis</i>	<i>Respiratory Syncytial Virus A</i>
<i>Pseudomonas aeruginosa</i>	<i>Respiratory Syncytial Virus B</i>
<i>Staphylococcus aureus</i>	<i>Rhinovirus</i>
<i>Staphylococcus epidermidis</i>	<i>SARS Coronavirus</i>
<i>Streptococcus pneumoniae</i>	<i>COVID-BEI GIV (SARS-CoV-2)</i>
<i>Streptococcus pyogenes</i>	<i>Candida albicans</i>
<i>Streptococcus salivarius</i>	<i>Pneumocystis jirovecii</i>

The complete performance characteristics of the TaqMan-SARS-Cov-2/Flu A/Flu B RT-PCR assay for 25 dry swabs per pool testing strategy is listed in **Table 19** below

Table 19: Performance characteristics of the TaqMan-SARS-Cov-2/Flu A/Flu B RT-PCR assay

Performance Metric		Description	Performance
Accuracy	PPA	TP/ (TP + FN)	269 / 270 = 99.62%; 95% CI [97.95% - 99.99%]
	PNA	TN/ (TN + FP)	90 / 90 = 100%; 95% CI [95.98% - 100%]
Precision	Repeatability	Within run replicate agreement	359 / 360 = 99.72%; 95% CI [98.46% - 99.99%]
	Reproducibility	Between run replicate agreement	119 / 120 = 99.16%; 95% CI [95.44% - 99.98%]
Analytical Sensitivity	Limit of Detection	Lowest concentration of virus detected with 95% probability of obtaining a correct result	1,500 copies/mL (ATCC VR-1986HK)
	Inclusivity	Amplification of all possible SARS-CoV-2/FluA/FluB genomes	100% homology to each of the sequences analyzed
	Cross-Reactivity	Amplification of non-target sequences	0% false positive and 0% false negative


		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report	
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989
Originating Department Analytics	Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 12 of 23

Validation Report and Data Review

All data and analysis were reviewed by laboratory director and validation report was compiled by the laboratory director listed below

Sameer Kalghatgi, Ph.D.
Sr. Director, Laboratory Operations
Infinity Biologix

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	


Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 13 of 23

Appendix: Detailed run data:

All raw data related to all the qPCR runs for obtaining the data in the appendix are stored securely in the ThermoFisher cloud platform under a common filename format – RUO-Gingko_Concentric Pooling_ExptDate. The cloud account is only accessible to authorized users with login credentials. Login information is changed every 6 months.


Day 1

IBX Sample	MS2	COVID N/S	IBX Result	Expected Call	COVID copies/mL
603299_1	15.1	32	Detected	Positive	1500
603299_10	20.87	34.67	Detected	Positive	1500
603299_2	18.76	31.85	Detected	Positive	1500
603299_3	19.56	21.72	Detected	Positive	1500
603299_4	18.6	33.93	Detected	Positive	1500
603299_5	16.29	32.95	Detected	Positive	1500
603299_6	17.37	32.28	Detected	Positive	1500
603299_7	16.98	32.14	Detected	Positive	1500
603299_8	20.34	33.95	Detected	Positive	1500
603299_9	19.86	34.22	Detected	Positive	1500
603300_1	19.07	34.14	Detected	Positive	1500
603300_10	19.31	28.98	Detected	Positive	1500
603300_2	19.37	33.85	Detected	Positive	1500
603300_3	19.23	34.38	Detected	Positive	1500
603300_4	19.78	35.39	Detected	Positive	1500
603300_5	19.68	34.34	Detected	Positive	1500
603300_6	19.84	34.59	Detected	Positive	1500
603300_7	19.95	34.95	Detected	Positive	1500
603300_8	19.37	32.46	Detected	Positive	1500
603300_9	19.19	29.54	Detected	Positive	1500
603302_1	19.48	29.12	Detected	Positive	1500
603302_10	19.24	29.62	Detected	Positive	1500
603302_2	19.38	27.88	Detected	Positive	1500
603302_3	19.54	28.51	Detected	Positive	1500
603302_4	18.83	29	Detected	Positive	1500
603302_5	18.57	28.69	Detected	Positive	1500
603302_6	18.11	31.56	Detected	Positive	1500
603302_7	18.44	30.86	Detected	Positive	1500
603302_8	18.7	31.17	Detected	Positive	1500
603302_9	19.49	31.19	Detected	Positive	1500
603297_1	19.93	31.86	Detected	Positive	1500
603297_10	18.9	30.51	Detected	Positive	1500

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	


Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 14 of 23

603297_2	19.6	30.85	Detected	Positive	1500
603297_3	19.79	31.52	Detected	Positive	1500
603297_4	19.8	30.75	Detected	Positive	1500
603297_5	20.56	32.97	Detected	Positive	1500
603297_6	20.07	31.74	Detected	Positive	1500
603297_7	20.1	31.82	Detected	Positive	1500
603297_8	20.1	30.98	Detected	Positive	1500
603297_9	18.5	30.43	Detected	Positive	1500
603301_1	19.23	25.93	Detected	Positive	1500
603301_10	20.17	30.3	Detected	Positive	1500
603301_2	19.23	24.86	Detected	Positive	1500
603301_3	20.09	28.99	Detected	Positive	1500
603301_4	19.17	29.3	Detected	Positive	1500
603301_5	18.82	25.67	Detected	Positive	1500
603301_6	18.47	29.45	Detected	Positive	1500
603301_7	20.01	30.3	Detected	Positive	1500
603301_8	19.67	31.66	Detected	Positive	1500
603301_9	20.22	32.4	Detected	Positive	1500
603303_1	18.6	29.08	Detected	Positive	1500
603303_10	18.88	31.75	Detected	Positive	1500
603303_2	19.31	22.43	Detected	Positive	1500
603303_3	18.33	26.95	Detected	Positive	1500
603303_4	19.45	25.41	Detected	Positive	1500
603303_5	19.91	24.45	Detected	Positive	1500
603303_6	19.68	26.45	Detected	Positive	1500
603303_7	19.11	24.12	Detected	Positive	1500
603303_8	19.4	26.51	Detected	Positive	1500
603303_9	19.49	28.58	Detected	Positive	1500
603304_1	18.35	30.52	Detected	Positive	1500
603304_10	20.88	26.75	Detected	Positive	1500
603304_2	19.41	31.86	Detected	Positive	1500
603304_3	17.87	32.2	Detected	Positive	1500
603304_4	18.22	28.79	Detected	Positive	1500
603304_5	19.88	22.67	Detected	Positive	1500
603304_6	21.22	27.89	Detected	Positive	1500
603304_7	20.38	30.42	Detected	Positive	1500
603304_8	20.81	25.95	Detected	Positive	1500
603304_9	21	28.12	Detected	Positive	1500
603298_1	20.14	20.39	Detected	Positive	1312500
603298_10	20.65	21.67	Detected	Positive	1312500

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 15 of 23

603298_2	21.22	22.11	Detected	Positive	1312500
603298_3	19.77	18.13	Detected	Positive	1312500
603298_4	19.39	20.12	Detected	Positive	1312500
603298_5	18.91	16.46	Detected	Positive	1312500
603298_6	19.47	21.19	Detected	Positive	1312500
603298_7	Undetermined	Undetermined	Inconclusive	Positive	1312500
603298_8	19.89	20.95	Detected	Positive	1312500
603298_9	19.6	20.12	Detected	Positive	1312500
603308_1	20.63	21.28	Detected	Positive	1312500
603308_10	20.41	20.59	Detected	Positive	1312500
603308_2	23.21	23.73	Detected	Positive	1312500
603308_3	21.05	21.78	Detected	Positive	1312500
603308_4	20.89	20.9	Detected	Positive	1312500
603308_5	20.46	21.05	Detected	Positive	1312500
603308_6	20.49	21.16	Detected	Positive	1312500
603308_7	20.95	21.27	Detected	Positive	1312500
603308_8	19.3	19.16	Detected	Positive	1312500
603308_9	20.56	21.07	Detected	Positive	1312500
603305_1	19.38	Undetermined	Not Detected	Negative	0
603305_10	18.67	Undetermined	Not Detected	Negative	0
603305_2	20.05	Undetermined	Not Detected	Negative	0
603305_3	19.22	Undetermined	Not Detected	Negative	0
603305_4	18.95	Undetermined	Not Detected	Negative	0
603305_5	19.15	Undetermined	Not Detected	Negative	0
603305_6	18.84	Undetermined	Not Detected	Negative	0
603305_7	18.87	Undetermined	Not Detected	Negative	0
603305_8	19.32	Undetermined	Not Detected	Negative	0
603305_9	19.24	Undetermined	Not Detected	Negative	0
603306_1	20.6	Undetermined	Not Detected	Negative	0
603306_10	19.81	Undetermined	Not Detected	Negative	0
603306_2	17.74	Undetermined	Not Detected	Negative	0
603306_3	14.55	Undetermined	Not Detected	Negative	0
603306_4	18.21	Undetermined	Not Detected	Negative	0
603306_5	16.95	Undetermined	Not Detected	Negative	0
603306_6	18.14	Undetermined	Not Detected	Negative	0
603306_7	18.1	Undetermined	Not Detected	Negative	0
603306_8	18.6	Undetermined	Not Detected	Negative	0
603306_9	19.89	Undetermined	Not Detected	Negative	0
603307_1	19.4	Undetermined	Not Detected	Negative	0
603307_10	20.08	Undetermined	Not Detected	Negative	0


		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 16 of 23

603307_2	19.47	Undetermined	Not Detected	Negative	0
603307_3	19.59	Undetermined	Not Detected	Negative	0
603307_4	19.04	Undetermined	Not Detected	Negative	0
603307_5	20.22	Undetermined	Not Detected	Negative	0
603307_6	19.54	Undetermined	Not Detected	Negative	0
603307_7	19.33	Undetermined	Not Detected	Negative	0
603307_8	19.68	Undetermined	Not Detected	Negative	0
603307_9	19.65	Undetermined	Not Detected	Negative	0


Day 2

IBX Sample	MS2	COVID N/S	IBX Result	Expected Call	COVID copies/mL
606490-1	22.36	33.24	Detected	Positive	1500
606490-10	22.11	31.30	Detected	Positive	1500
606490-2	22.21	33.17	Detected	Positive	1500
606490-3	22.00	33.49	Detected	Positive	1500
606490-4	22.20	33.67	Detected	Positive	1500
606490-5	22.21	33.84	Detected	Positive	1500
606490-6	22.29	33.97	Detected	Positive	1500
606490-7	22.33	33.89	Detected	Positive	1500
606490-8	22.06	33.59	Detected	Positive	1500
606490-9	22.43	33.94	Detected	Positive	1500
606491-1	21.84	32.38	Detected	Positive	1500
606491-10	21.97	32.90	Detected	Positive	1500
606491-2	22.00	33.67	Detected	Positive	1500
606491-3	22.27	33.06	Detected	Positive	1500
606491-4	21.98	33.06	Detected	Positive	1500
606491-5	21.85	33.93	Detected	Positive	1500
606491-6	21.47	34.13	Detected	Positive	1500
606491-7	22.30	33.12	Detected	Positive	1500
606491-8	22.00	33.58	Detected	Positive	1500
606491-9	21.97	33.34	Detected	Positive	1500
606494-1	22.29	33.98	Detected	Positive	1500
606494-10	22.03	34.39	Detected	Positive	1500
606494-2	22.53	33.77	Detected	Positive	1500
606494-3	22.28	34.00	Detected	Positive	1500
606494-4	22.25	33.91	Detected	Positive	1500
606494-5	22.37	34.00	Detected	Positive	1500
606494-6	22.28	33.50	Detected	Positive	1500

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	


Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 17 of 23

606494-7	22.29	34.00	Detected	Positive	1500
606494-8	22.29	34.09	Detected	Positive	1500
606494-9	22.24	33.82	Detected	Positive	1500
606488-1	21.80	30.57	Detected	Positive	1500
606488-10	21.83	30.12	Detected	Positive	1500
606488-2	22.27	30.83	Detected	Positive	1500
606488-3	22.09	31.26	Detected	Positive	1500
606488-4	22.07	30.75	Detected	Positive	1500
606488-5	22.06	30.53	Detected	Positive	1500
606488-6	22.39	31.50	Detected	Positive	1500
606488-7	22.11	30.42	Detected	Positive	1500
606488-8	22.16	30.49	Detected	Positive	1500
606488-9	21.98	30.54	Detected	Positive	1500
606495-1	22.44	30.55	Detected	Positive	1500
606495-10	22.75	31.33	Detected	Positive	1500
606495-2	22.46	30.70	Detected	Positive	1500
606495-3	22.21	30.28	Detected	Positive	1500
606495-4	22.08	30.53	Detected	Positive	1500
606495-5	22.20	30.49	Detected	Positive	1500
606495-6	21.86	30.39	Detected	Positive	1500
606495-7	21.95	30.14	Detected	Positive	1500
606495-8	22.06	30.41	Detected	Positive	1500
606495-9	22.10	30.44	Detected	Positive	1500
606496-1	22.16	31.12	Detected	Positive	1500
606496-10	22.04	30.71	Detected	Positive	1500
606496-2	23.07	33.46	Detected	Positive	1500
606496-3	22.58	31.04	Detected	Positive	1500
606496-4	23.47	35.13	Detected	Positive	1500
606496-5	22.30	32.20	Detected	Positive	1500
606496-6	22.37	32.10	Detected	Positive	1500
606496-7	22.06	30.81	Detected	Positive	1500
606496-8	22.24	30.92	Detected	Positive	1500
606496-9	22.27	31.15	Detected	Positive	1500
606497-1	22.06	30.65	Detected	Positive	1500
606497-10	21.80	30.32	Detected	Positive	1500
606497-2	21.90	30.44	Detected	Positive	1500
606497-3	22.06	30.26	Detected	Positive	1500
606497-4	21.98	30.52	Detected	Positive	1500
606497-5	22.02	31.07	Detected	Positive	1500

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 18 of 23

606497-6	22.35	30.70	Detected	Positive	1500
606497-7	21.84	30.43	Detected	Positive	1500
606497-8	22.23	30.58	Detected	Positive	1500
606497-9	22.04	31.20	Detected	Positive	1500
606489-1	21.88	20.88	Detected	Positive	1312500
606489-10	21.91	20.71	Detected	Positive	1312500
606489-2	22.18	21.50	Detected	Positive	1312500
606489-3	21.71	20.73	Detected	Positive	1312500
606489-4	22.91	21.23	Detected	Positive	1312500
606489-5	21.67	20.88	Detected	Positive	1312500
606489-6	21.91	20.77	Detected	Positive	1312500
606489-7	21.80	20.78	Detected	Positive	1312500
606489-8	22.60	21.73	Detected	Positive	1312500
606489-9	21.87	20.78	Detected	Positive	1312500
606498-1	22.34	20.76	Detected	Positive	1312500
606498-10	22.44	20.86	Detected	Positive	1312500
606498-2	22.38	20.89	Detected	Positive	1312500
606498-3	22.33	20.82	Detected	Positive	1312500
606498-4	22.95	20.85	Detected	Positive	1312500
606498-5	22.15	20.73	Detected	Positive	1312500
606498-6	22.25	20.71	Detected	Positive	1312500
606498-7	22.30	20.67	Detected	Positive	1312500
606498-8	22.41	20.83	Detected	Positive	1312500
606498-9	22.43	20.97	Detected	Positive	1312500
606492-1	22.19	Undetermined	Not Detected	Negative	0
606492-10	21.84	Undetermined	Not Detected	Negative	0
606492-2	21.63	Undetermined	Not Detected	Negative	0
606492-3	22.50	Undetermined	Not Detected	Negative	0
606492-4	21.98	Undetermined	Not Detected	Negative	0
606492-5	21.62	Undetermined	Not Detected	Negative	0
606492-6	21.50	Undetermined	Not Detected	Negative	0
606492-7	21.63	Undetermined	Not Detected	Negative	0
606492-8	21.61	Undetermined	Not Detected	Negative	0
606492-9	22.15	Undetermined	Not Detected	Negative	0
606493-1	21.80	Undetermined	Not Detected	Negative	0
606493-10	22.65	Undetermined	Not Detected	Negative	0
606493-2	21.82	Undetermined	Not Detected	Negative	0
606493-3	21.32	Undetermined	Not Detected	Negative	0
606493-4	22.05	Undetermined	Not Detected	Negative	0


		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 19 of 23

606493-5	21.62	Undetermined	Not Detected	Negative	0
606493-6	21.74	Undetermined	Not Detected	Negative	0
606493-7	22.00	Undetermined	Not Detected	Negative	0
606493-8	21.87	Undetermined	Not Detected	Negative	0
606493-9	22.09	Undetermined	Not Detected	Negative	0
606499-1	22.16	Undetermined	Not Detected	Negative	0
606499-10	22.13	Undetermined	Not Detected	Negative	0
606499-2	22.06	Undetermined	Not Detected	Negative	0
606499-3	22.10	Undetermined	Not Detected	Negative	0
606499-4	21.88	Undetermined	Not Detected	Negative	0
606499-5	22.07	Undetermined	Not Detected	Negative	0
606499-6	22.27	Undetermined	Not Detected	Negative	0
606499-7	22.32	Undetermined	Not Detected	Negative	0
606499-8	22.06	Undetermined	Not Detected	Negative	0
606499-9	22.46	Undetermined	Not Detected	Negative	0


Day 3

IBX Sample	MS2	COVID N/S	IBX Result	Expected Call	COVID copies/mL
608813-1	23.0	34.5	Detected	Positive	1500
608813-10	22.5	29.8	Detected	Positive	1500
608813-2	22.6	33.4	Detected	Positive	1500
608813-3	22.8	34.5	Detected	Positive	1500
608813-4	22.3	33.1	Detected	Positive	1500
608813-5	22.4	32.5	Detected	Positive	1500
608813-6	22.2	32.7	Detected	Positive	1500
608813-7	22.6	33.1	Detected	Positive	1500
608813-8	22.3	33.6	Detected	Positive	1500
608813-9	22.3	34.0	Detected	Positive	1500
608815-1	22.0	33.7	Detected	Positive	1500
608815-10	22.7	29.7	Detected	Positive	1500
608815-2	22.3	35.3	Detected	Positive	1500
608815-3	22.9	35.1	Detected	Positive	1500
608815-4	22.7	32.0	Detected	Positive	1500
608815-5	22.1	32.9	Detected	Positive	1500
608815-6	22.5	34.3	Detected	Positive	1500
608815-7	22.0	33.7	Detected	Positive	1500
608815-8	22.8	34.4	Detected	Positive	1500
608815-9	22.7	34.3	Detected	Positive	1500
608821-1	23.8	32.6	Detected	Positive	1500

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	


Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 20 of 23

608821-10	22.9	34.0	Detected	Positive	1500
608821-2	22.2	23.5	Detected	Positive	1500
608821-3	23.5	28.6	Detected	Positive	1500
608821-4	22.9	29.7	Detected	Positive	1500
608821-5	22.6	31.5	Detected	Positive	1500
608821-6	22.1	32.4	Detected	Positive	1500
608821-7	23.1	34.2	Detected	Positive	1500
608821-8	22.8	32.9	Detected	Positive	1500
608821-9	22.1	32.6	Detected	Positive	1500
608811-1	22.6	33.7	Detected	Positive	1500
608811-10	22.4	33.1	Detected	Positive	1500
608811-2	22.3	33.2	Detected	Positive	1500
608811-3	21.9	26.9	Detected	Positive	1500
608811-4	21.7	30.3	Detected	Positive	1500
608811-5	22.4	32.6	Detected	Positive	1500
608811-6	22.2	32.0	Detected	Positive	1500
608811-7	22.3	32.8	Detected	Positive	1500
608811-8	22.8	32.6	Detected	Positive	1500
608811-9	22.8	34.2	Detected	Positive	1500
608818-1	22.4	26.5	Detected	Positive	1500
608818-10	22.8	33.5	Detected	Positive	1500
608818-2	22.6	27.6	Detected	Positive	1500
608818-3	22.5	29.4	Detected	Positive	1500
608818-4	23.2	32.4	Detected	Positive	1500
608818-5	22.6	29.8	Detected	Positive	1500
608818-6	22.5	28.7	Detected	Positive	1500
608818-7	22.8	34.0	Detected	Positive	1500
608818-8	22.4	33.7	Detected	Positive	1500
608818-9	22.8	33.2	Detected	Positive	1500
608819-1	23.7	34.8	Detected	Positive	1500
608819-10	22.6	31.8	Detected	Positive	1500
608819-2	23.3	34.4	Detected	Positive	1500
608819-3	23.2	35.0	Detected	Positive	1500
608819-4	22.2	32.6	Detected	Positive	1500
608819-5	23.2	34.1	Detected	Positive	1500
608819-6	24.3	38.0	Detected	Positive	1500
608819-7	22.0	32.9	Detected	Positive	1500
608819-8	22.0	32.3	Detected	Positive	1500
608819-9	22.3	31.7	Detected	Positive	1500

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	


Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 21 of 23

608820-1	22.8	30.7	Detected	Positive	1500
608820-10	22.1	29.8	Detected	Positive	1500
608820-2	22.6	31.3	Detected	Positive	1500
608820-3	23.5	33.2	Detected	Positive	1500
608820-4	22.7	32.6	Detected	Positive	1500
608820-5	22.9	32.6	Detected	Positive	1500
608820-6	23.1	24.9	Detected	Positive	1500
608820-7	22.8	32.7	Detected	Positive	1500
608820-8	23.3	32.4	Detected	Positive	1500
608820-9	23.2	30.6	Detected	Positive	1500
608814-1	23.4	23.6	Detected	Positive	1312500
608814-10	22.9	21.8	Detected	Positive	1312500
608814-2	22.9	21.6	Detected	Positive	1312500
608814-3	22.4	21.5	Detected	Positive	1312500
608814-4	22.7	22.9	Detected	Positive	1312500
608814-5	22.8	21.7	Detected	Positive	1312500
608814-6	23.1	21.6	Detected	Positive	1312500
608814-7	22.8	22.8	Detected	Positive	1312500
608814-8	25.3	26.8	Detected	Positive	1312500
608814-9	22.9	22.7	Detected	Positive	1312500
608822-1	22.6	22.1	Detected	Positive	1312500
608822-10	23.4	22.6	Detected	Positive	1312500
608822-2	22.8	22.7	Detected	Positive	1312500
608822-3	22.5	22.4	Detected	Positive	1312500
608822-4	22.3	21.8	Detected	Positive	1312500
608822-5	22.4	21.6	Detected	Positive	1312500
608822-6	22.4	21.4	Detected	Positive	1312500
608822-7	23.2	21.6	Detected	Positive	1312500
608822-8	22.0	20.7	Detected	Positive	1312500
608822-9	23.0	22.2	Detected	Positive	1312500
608816-1	23.5	Undetermined	Not Detected	Negative	0
608816-10	23.3	Undetermined	Not Detected	Negative	0
608816-2	23.6	Undetermined	Not Detected	Negative	0
608816-3	23.7	Undetermined	Not Detected	Negative	0
608816-4	22.8	Undetermined	Not Detected	Negative	0
608816-5	23.4	Undetermined	Not Detected	Negative	0
608816-6	23.5	Undetermined	Not Detected	Negative	0
608816-7	23.2	Undetermined	Not Detected	Negative	0
608816-8	22.8	Undetermined	Not Detected	Negative	0

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report			
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989		
Originating Department Analytics		Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey	

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 22 of 23

608816-9	22.7	Undetermined	Not Detected	Negative	0
608812-1	22.6	Undetermined	Not Detected	Negative	0
608812-10	22.4	Undetermined	Not Detected	Negative	0
608812-2	22.7	Undetermined	Not Detected	Negative	0
608812-3	22.4	Undetermined	Not Detected	Negative	0
608812-4	22.6	Undetermined	Not Detected	Negative	0
608812-5	22.2	Undetermined	Not Detected	Negative	0
608812-6	22.7	Undetermined	Not Detected	Negative	0
608812-7	22.8	Undetermined	Not Detected	Negative	0
608812-8	22.9	Undetermined	Not Detected	Negative	0
608812-9	22.6	Undetermined	Not Detected	Negative	0
608817-1	22.4	Undetermined	Not Detected	Negative	0
608817-10	22.3	Undetermined	Not Detected	Negative	0
608817-2	22.6	Undetermined	Not Detected	Negative	0
608817-3	22.1	Undetermined	Not Detected	Negative	0
608817-4	22.7	Undetermined	Not Detected	Negative	0
608817-5	21.5	Undetermined	Not Detected	Negative	0
608817-6	25.0	Undetermined	Not Detected	Negative	0
608817-7	22.5	Undetermined	Not Detected	Negative	0
608817-8	23.7	Undetermined	Not Detected	Negative	0
608817-9	22.0	Undetermined	Not Detected	Negative	0

		Category: 3003 Validation Records Title: IBX Gingko Swab Pooling Surveillance Assay Validation Report	
Version 03	State Effective	Effective Date 03-MAR-2021	Document ID 441989
Originating Department Analytics	Affected Department(s) Analytics		Affected Site(s) IBX-New Jersey

Printed by Sameer.Kalghatgi@Infinity-biologix.com from app.zenqms.com on 03-Mar-2021 at 7:13:00 PM UTC • Page 23 of 23

REVISION HISTORY

Version 01 Effective on 19-Feb-2021
New Document

Version 02 Effective on 02-Mar-2021
Minor update to language

Version 03 Effective on 03-Mar-2021
Document category changed

DOCUMENT ELECTRONIC SIGNATURES

DOCUMENT APPROVAL WORKFLOW

Author Approval

Sameer Kalghatgi
Sr. Director Lab Operations
Sameer.Kalghatgi@Infinity-biologix.com

I am the author of this document.
Signed 4:24:35 PM UTC 03-Mar-2021

Required Workflow Steps for this Category

Christian Bixby
Senior Director, Analytical and Clinical Services
Christian.Bixby@infinity-biologix.com

Infinity Biologix / Approver 1
I have reviewed and approve this document.
Signed 4:36:57 PM UTC 03-Mar-2021

Shareef Nahas
VP, Global Lab Operations
shareef.nahas@infinity-biologix.com

Infinity Biologix / Approver 2
I have reviewed and approve this document.
Signed 6:33:47 PM UTC 03-Mar-2021

Mayur Patel
Sr. Director, Quality Assurance
Mayur.Patel@infinity-biologix.com

Infinity Biologix / Quality Approval
I have reviewed and approve this document.
Signed 6:40:23 PM UTC 03-Mar-2021